

Proteomics

What is proteomics?

- Studying the protein content of a particular organism, organ or cell type
- We have to resolve all the proteins in a particular system
- This can number over 40 000 proteins
- Proteomics is part of functional genomics

Why use Proteomics ?

- Allows us to compare control and experimental samples to look for differences in the protein profile
- Looks directly at the protein level rather than trying to make inferences from genetic information
- We can eventually use the information gathered to isolate the genes involved


Genomics can be limited.

- From a protein function aspect there maybe post-translational modifications. These modifications may not be apparent from gene sequence
 - Phosphorylation
 - Glycosylation
 - Germylation (C20) Membrane assoc.
 - Palmitoylation (C16) Membrane assoc.

Proteomics

- Divided up into two distinct areas
 - 2-D gel electrophoresis for the characterisation of protein profiles
 - Mass spectroscopy for the accurate identification of isolated proteins, including post-translational modifications.

Proteins are complicated macro-molecules



We can be reductionist and look at amino acids....





table 5-1

| Properties and Conventions Associated with the Standard Amino Acids | | | | | | | | | |
|---|-------------------|---|-------|-------------------|--|---------------------|-------|-------------------------------|---|
| Amino acid | Abbreviated names | | M_r | pK_a values | | | pI | Hydropathy index [*] | Occurrence in proteins (%) [†] |
| | | | | pK_1 (—COOH) | pK_2 (—NH ₃ ⁺) | pK_R (R group) | | | |
| Nonpolar, aliphatic R groups | | | | | | | | | |
| Glycine | Gly | G | 75 | 2.34 | 9.60 | | 5.97 | -0.4 | 7.2 |
| Alanine | Ala | A | 89 | 2.34 | 9.69 | | 6.01 | 1.8 | 7.8 |
| Valine | Val | V | 117 | 2.32 | 9.62 | | 5.97 | 4.2 | 6.6 |
| Leucine | Leu | L | 131 | 2.36 | 9.60 | | 5.98 | 3.8 | 9.1 |
| Isoleucine | Ile | I | 131 | 2.36 | 9.68 | | 6.02 | 4.5 | 5.3 |
| Methionine | Met | M | 149 | 2.28 | 9.21 | | 5.74 | 1.9 | 2.3 |
| Aromatic R groups | | | | | | | | | |
| Phenylalanine | Phe | F | 165 | 1.83 | 9.13 | | 5.48 | 2.8 | 3.9 |
| Tyrosine | Tyr | Y | 181 | 2.20 | 9.11 | 10.07 | 5.66 | -1.3 | 3.2 |
| Tryptophan | Trp | W | 204 | 2.38 | 9.39 | | 5.89 | -0.9 | 1.4 |
| Polar, uncharged R groups | | | | | | | | | |
| Serine | Ser | S | 105 | 2.21 | 9.15 | | 5.68 | -0.8 | 6.8 |
| Proline | Pro | P | 115 | 1.99 | 10.96 | | 6.48 | 1.6 | 5.2 |
| Threonine | Thr | T | 119 | 2.11 | 9.62 | | 5.87 | -0.7 | 5.9 |
| Cysteine | Cys | C | 121 | 1.96 | 10.28 | 8.18 | 5.07 | 2.5 | 1.9 |
| Asparagine | Asn | N | 132 | 2.02 | 8.80 | | 5.41 | -3.5 | 4.3 |
| Glutamine | Gln | Q | 146 | 2.17 | 9.13 | | 5.65 | -3.5 | 4.2 |
| Positively charged R groups | | | | | | | | | |
| Lysine | Lys | K | 146 | 2.18 | 8.95 | 10.53 | 9.74 | -3.9 | 5.9 |
| Histidine | His | H | 155 | 1.82 | 9.17 | 6.00 | 7.59 | -3.2 | 2.3 |
| Arginine | Arg | R | 174 | 2.17 | 9.04 | 12.48 | 10.76 | -4.5 | 5.1 |
| Negatively charged R groups | | | | | | | | | |
| Aspartate | Asp | D | 133 | 1.88 | 9.60 | 3.65 | 2.77 | -3.5 | 5.3 |
| Glutamate | Glu | E | 147 | 2.19 | 9.67 | 4.25 | 3.22 | -3.5 | 6.3 |


*A scale combining hydrophobicity and hydrophilicity of R groups; it can be used to measure the tendency of an amino acid to seek an aqueous environment (— values) or a hydrophobic environment (+ values). See Chapter 12. From Kyte, J. & Doolittle, R.F. (1982) *J. Mol. Biol.* **157**, 105–132.

[†]Average occurrence in over 1150 proteins. From Doolittle, R.F. (1989) Redundancies in protein sequences. In *Prediction of Protein Structure and the Principles of Protein Conformation* (Fasman, G.D., ed) Plenum Press, NY, pp. 599–623.

In effect all proteins are zwitterions




Nonionic
form



Zwitterionic
form

Proteins are sequential molecules with defined ends



Proteins vary greatly in size

table 5–2

Molecular Data on Some Proteins

| | Molecular weight | Number of residues | Number of polypeptide chains |
|---|------------------|--------------------|------------------------------|
| Cytochrome <i>c</i> (human) | 13,000 | 104 | 1 |
| Ribonuclease A (bovine pancreas) | 13,700 | 124 | 1 |
| Lysozyme (egg white) | 13,930 | 129 | 1 |
| Myoglobin (equine heart) | 16,890 | 153 | 1 |
| Chymotrypsin (bovine pancreas) | 21,600 | 241 | 3 |
| Chymotrypsinogen (bovine) | 22,000 | 245 | 1 |
| Hemoglobin (human) | 64,500 | 574 | 4 |
| Serum albumin (human) | 68,500 | 609 | 1 |
| Hexokinase (yeast) | 102,000 | 972 | 2 |
| RNA polymerase (<i>E. coli</i>) | 450,000 | 4,158 | 5 |
| Apolipoprotein B (human) | 513,000 | 4,536 | 1 |
| Glutamine synthetase (<i>E. coli</i>) | 619,000 | 5,628 | 12 |
| Titin (human) | 2,993,000 | 26,926 | 1 |

Amino acid composition defines the physical nature of a protein

- Each amino acid has a particular pI
- These pI's combine to give a protein an over all pI.
- As proteins are zwitterions when they are maintained at their pI they have no charge.

table 5-3

| Amino Acid Composition of Two Proteins* | | |
|---|--|-------------------------|
| Amino acid | Number of residues per molecule of protein | |
| | Bovine cytochrome c | Bovine chymotrypsinogen |
| Ala | 6 | 22 |
| Arg | 2 | 4 |
| Asn | 5 | 15 |
| Asp | 3 | 8 |
| Cys | 2 | 10 |
| Gln | 3 | 10 |
| Glu | 9 | 5 |
| Gly | 14 | 23 |
| His | 3 | 2 |
| Ile | 6 | 10 |
| Leu | 6 | 19 |
| Lys | 18 | 14 |
| Met | 2 | 2 |
| Phe | 4 | 6 |
| Pro | 4 | 9 |
| Ser | 1 | 28 |
| Thr | 8 | 23 |
| Trp | 1 | 8 |
| Tyr | 4 | 4 |
| Val | 3 | 23 |
| Total | 104 | 245 |

*Note that standard procedures for the acid hydrolysis of proteins convert Asn and Gln to Asp and Glu, respectively. In addition, Trp is destroyed. Special procedures must be employed to determine the amounts of these amino acids.

Protein separation for proteomic investigations

1st dimension

Separation by charge


Proteins vary in isoelectric point

table 5–6

| The Isoelectric Points of Some Proteins | |
|---|------|
| Protein | pI |
| Pepsin | ~1.0 |
| Egg albumin | 4.6 |
| Serum albumin | 4.9 |
| Urease | 5.0 |
| β -Lactoglobulin | 5.2 |
| Hemoglobin | 6.8 |
| Myoglobin | 7.0 |
| Chymotrypsinogen | 9.5 |
| Cytochrome <i>c</i> | 10.7 |
| Lysozyme | 11.0 |

Isoelectric focusing.

Separation by charge (pI)



Protein separation for proteomic investigations


2nd dimension

Separation by size

Sodium dodecyl sulphate polyacrylamide electrophoresis (SDS-PAGE)

- System relies on two chemicals
 - Acrylamide
 - Forms a gel-like matrix of known composition
 - The “denser” a gel is made, the slower proteins will move.
 - Sodium dodecyl sulphate
 - A detergent which binds to the protein giving all proteins a gross negative charge.
 - All proteins now migrate towards a +ve charge.

SDS gives all proteins a net negative charge




Sodium dodecyl sulfate
(SDS)

Acrylamide reaction


Acrylamide gels are formed by co-polymerisation of acrylamide monomer , $\text{CH}_2=\text{CH}-\text{CO}-\text{NH}_2$, and a crosslinking co-monomer, N,N'-methylenebisacrylamide, $\text{CH}_2=\text{CH}-\text{CO}-\text{NH}-\text{CH}_2-\text{NH}-\text{CO}-\text{CH}=\text{CH}_2$, bisacrylanide

The mechanism of gel formation is vinyl addition polymerisation and is catalysed by free radical-generation system composed of ammonium persulphate (initiator) and accelerator tetramethylethylenediamine (TEMED).

Components of acrylamide




Acrylamide




N,N'-methylene
bisacrylamide


Acrylamide polymerisation



When run on a gel all proteins will move dependant on molecular weight alone




Example of protein separation on PA gel




(b)

Sizing of proteins




(a)




(b)

2-D gel electrophoresis




Visualisation of proteins


- Coomassie stain
 - Quick (1 hour) and cheap (20p/gel)
 - Not very sensitive (30ng/spot)
- Fluorescent stains (SYPRO Orange)
 - Slower (3 hours)
 - Sensitive (5ng/spot)
- Silver stain
 - Slower (2 hours)
 - Very sensitive (0.5ng/spot)





From Gel Electrophoresis - Principles and Basics



Cy3




Cy5




overlap

A range of IEF gels can be used




- We now have a pH range from 2 to 10 over three gels.

Gel arrays can increase separation




Soy storage proteins separated by the ReadyStrip narrow overlapping IPG strips, courtesy of Dr. Jim Malone, Monsanto, Inc.

We can also carry out differential extractions to increase resolution.



We can compare samples from control and experimental systems



We are now able to isolate proteins whose expression patterns are altered by a particular treatment.

What to do with these spots?

- Image analysis and associated software can allow differences to be noted.
- This information can then be fed to a robot spot cutter.
- The stained protein spots are isolated and placed in 96 well plates.
- The spots are digested with trypsin.
- Peptide mass maps can allow the identification of proteins.

We can fragment proteins at specific amino acid residues

table 5-7


The Specificity of Some Common Methods for Fragmenting Polypeptide Chains

| Treatment* | Cleavage points† |
|---|-------------------|
| Trypsin | Lys, Arg (C) |
| <i>Submaxillaris</i> protease | Arg (C) |
| Chymotrypsin | Phe, Trp, Tyr (C) |
| <i>Staphylococcus aureus</i> V8 protease | Asp, Glu (C) |
| Asp-N-protease | Asp, Glu (N) |
| Pepsin | Phe, Trp, Tyr (N) |
| Endoproteinase Lys C | Lys (C) |
| Cyanogen bromide | Met (C) |

*All except cyanogen bromide are proteases. All are available from commercial sources.

†Residues furnishing the primary recognition point for the protease or reagent; peptide bond cleavage occurs on either the carbonyl (C) or the amino (N) side of the indicated amino acid residues.

Cleavage of proteins can allow sequence determination




Once the amino acid sequence is known the gene can be isolated

| | |
|----------------------------------|-------------------------|
| Amino acid sequence (protein) | Gln–Tyr–Pro–Thr–Ile–Trp |
| | └──┴──┴──┴──┴──┴──┴──┘ |
| DNA sequence (gene) | CAGTATCCTACGATTTGG |


How Does MALDI-TOF Work?

- The sample is mixed with the matrix and an aliquot is deposited onto the target plate.
- The target plate spot where the sample is located is irradiated by a pulse of light from the laser.
- Ions are desorbed into the gas phase.
- Ions are accelerated to a high kinetic energy by the plates of the ion source. These ions enter the flight tube. The unit of time that it takes the ion to move from the source to the detector is called the TIME-OF-FLIGHT.
- Ions strike the detector which converts them to electrical signal.

Mass spectroscopy





Varying the magnetic field allows different ions to pass through to the detector.




(b)

How does electrospray MS-MS work?

- A mixture of peptides (in our case the tryptic digested spot) are sprayed through a fine needle (1 μ m) into the mass spec.
- Individual peptides are isolated during the first step and then fragmented during the second step.
- Fragments are obtained from the N- or C-terminal and are designated 'b' or 'y' ions.
- It is then a case of using software to look at ion profiles and reconstruct the protein sequence.



(a)




(b)

How do these two systems work together?

- We can now isolate proteins from a 2-D gel and identify the protein using mass spectroscopy.
- This type of process can be carried out rapidly allowing the development of high through-put systems.
- We can identify up to 200 spots a day from peptide mass maps and 20 from MS-MS.


Orbitrap

Ion trap mass analyzer consisting of an outer barrel-like electrode and a coaxial inner spindle-like electrode that traps ions in an orbital motion around the spindle.




Orbitrap


The image of trapped ions is detected and converted to a mass spectrum by first using the Fourier transform of time domain of the harmonic to create a frequency signal which is converted to mass



Top-down and bottom-up sequencing for identifying a specific protein or proteins of interest



Summary slide



https://www.matrixscience.com/cgi/search_form.pl?FORMVER=2&SEARCH=PMF

https://www.matrixscience.com/training/pmf1_q.html